


EQUIPMENT: OUTER BOOTS, WATER, OUTER GLOVES, INNER GLOVES, MASK, GOGGLES, SUIT



REMOVE ALL JEWELRY AND PERSONAL EFFECTS



CAREFULLY INSPECT THE SEAMS OF THE SUIT FOR TEARS OR BROKEN SEALS



STEP INTO SUIT AND PULL UP TO WAIST LEVEL

UNIVERSITY OF NEBRASKA MEDICAL CENTER^{**}

THEN APPLY OUTER BOOTS



TAPE BOOTS ALONG UPPER RIM SEALING THE BOOT TO THE SUIT



CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING







APPLY THE INNER GLOVES



APPLY THE UPPER PORTION OF THE SUIT.

UNIVERSITY OF NEBRASKA MEDICAL CENTER" HEROES

MAKE SURE SLEEVES OVERLAP INNER GLOVES





APPLY THE OUTER GLOVES



UNIVERSITY OF NEBRASKA MEDICAL CENTER" HEROES

TAPE GLOVE ALONG UPPER RIM SEALING THE GLOVE TO THE SUIT



CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING



REPEAT TAPING PROCESS FOR OTHER GLOVE



REMOVE TAPE FROM FRONT ZIPPER FLAP OF SUIT



SEAL ZIPPER FLAP





PULL THE HOOD OF SUIT UP AND MAKE FINAL CHECK OF TAPE SEALS

REMOVING SUIT SUPPORT PPE

- Begin by cleaning any areas of the suit that may have been chemically contaminated if necessary.
- Carefully remove all tape from gloves and boots and unzip suit.
- Remove outer gloves being careful not to contaminate inner gloves.
- Remove suit only touching inner surfaces. Try to keep outer surfaces turned in on themselves.
- Step into a clean (cold zone) area.
- Once suit is removed, inner gloves can be removed and hands sanitized if desired.
- Remove respiratory protection last, then proceed to hand washing and shower.

NOTE: If any areas of your body were contaminated during the doffing process, decontaminate them well in the shower and **notify the decon leader or appropriate personnel**.



UNIVERSITY OF NEBRASH MEDICAL CENTER" HEROES