

Type B - Donning





CAREFULLY
INSPECT THE
SEAMS OF THE
SUIT FOR TEARS
OR BROKEN
SEALS



DONNING PARTNER TAKES AND RECORDS VITAL SIGNS



REMOVE ALL JEWELRY AND PERSONAL EFFECTS



DRINK 1-2 GLASSES OF A NON-CAFFEINATED BEVERAGE



APPLY BOOT LINERS TO BOTH FEET AND TAPE TO OUTER SURFACE OF PANTS





Type B - Donning




CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING



STEP INTO SUIT AND PULL UP TO WAIST LEVEL







TAPE BOOTS ALONG UPPER RIM SEALNG THEM TO THE SUIT



CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING



APPLY INNER GLOVES





Type B - Donning



STAND AND APPLY THE UPPER PORTION OF THE SUIT



APPLY THE
OUTER GLOVES.

CAREFULLY
TUCK THE CUFF
OF THE SUIT
INSIDE



TAPE GLOVE ALONG UPPER RIM SEALING THE GLOVE TO THE SUIT



CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING




THE ZIPPER IS SEALED HALF WAY UP THE SUIT





Type B - Donning









PUT ON THE MASK



TIGHTEN
EACH STRAP
SYSTEMATICALLY, THEN
CHECK THE
SEAL

(FITTESTING SHOULD BE PERFORMED ANUALLY)



PUT THE HOOD UP AND CLOSE THE REST OF THE ZIPPER





Type B - Donning







CLICK REGULATOR **INTO PLACE** ON MASK. AIR FLOW SHOULD **START**



THE SUPPORT TEAM SHOULD BE AWARE OF HAND SIGNALS USED WHILE IN THE SUIT



PATTING THE HEAD **INDICATES NEEDING OUT** OF THE SUIT



PUTTING HANDS AT THE **NECK SIGNIFIES INABILITY TO BREATHE**

