



CHEMICAL PERSONAL PROTECTIVE EQUIPMENT

Type B - Donning



EQUIPMENT: AIR SUPPLY HOSE, REGULATOR, MASK, BELT, TAPE, SUIT, INNER BOOTS, OUTER BOOTS, INNER GLOVES, OUTER GLOVES



CAREFULLY INSPECT THE SEAMS OF THE SUIT FOR TEARS OR BROKEN SEALS



DONNING PARTNER TAKES AND RECORDS VITAL SIGNS



REMOVE ALL JEWELRY AND PERSONAL EFFECTS



DRINK 1-2 GLASSES OF A NON-CAFFEINATED BEVERAGE



APPLY BOOT LINERS TO BOTH FEET AND TAPE TO OUTER SURFACE OF PANTS





CHEMICAL PERSONAL PROTECTIVE EQUIPMENT

Type B - Donning

07



CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING

08



STEP INTO SUIT AND PULL UP TO WAIST LEVEL

09



APPLY OUTER BOOTS OVER THE SUIT. THE RUBBER GRIPS SHOULD BE ALONG THE BOTTOM TO IMPROVE TRACTION

10



TAPE BOOTS ALONG UPPER RIM SEALING THEM TO THE SUIT

11



CREATE A TAB ON THE TAPE'S END TO AID IN LATER DOFFING

12



APPLY INNER GLOVES





CHEMICAL PERSONAL PROTECTIVE EQUIPMENT

Type B - Donning

13



STAND AND
APPLY THE
UPPER PORTION
OF THE SUIT

14



APPLY THE
OUTER GLOVES.

CAREFULLY
TUCK THE CUFF
OF THE SUIT
INSIDE

15



TAPE GLOVE
ALONG UPPER
RIM SEALING
THE GLOVE TO
THE SUIT

16



CREATE A TAB
ON THE TAPE'S
END TO AID IN
LATER DOFFING

17



THE ZIPPER IS
SEALED HALF
WAY UP THE
SUIT



CHEMICAL PERSONAL PROTECTIVE EQUIPMENT

Type B - Donning



18

PUT ON THE BELT



19

BELT MOUNTED REGULATOR CLIP



20

BELT MOUNTED AIR SUPPLY HOSE CLIP



21

PUT ON THE MASK



22

TIGHTEN EACH STRAP SYSTEMATICALLY, THEN CHECK THE SEAL

(FIT TESTING SHOULD BE PERFORMED ANUALLY)



23

PUT THE HOOD UP AND CLOSE THE REST OF THE ZIPPER





CHEMICAL PERSONAL PROTECTIVE EQUIPMENT

Type B - Donning



24

TAPE NECK REGION FOR ADDED SPLASH PROTECTION



25

CONNECT AIR SUPPLY HOSE TO REGULATOR HOSE AT THE BELT. CHECK REGULATOR VALVES



26

CLICK REGULATOR INTO PLACE ON MASK. AIR FLOW SHOULD START



IMPORTANT

THE SUPPORT TEAM SHOULD BE AWARE OF HAND SIGNALS USED WHILE IN THE SUIT



!

PATting THE HEAD INDICATES NEEDING OUT OF THE SUIT



!

PUTTING HANDS AT THE NECK SIGNIFIES INABILITY TO BREATHE

